Medefil, Inc.
0.9% SODIUM CHLORIDE INJECTION, USP
"NORMAL SALINE I.V. FLUSH SYRINGE"

DESCRIPTION: 0.9% Sodium Chloride Injection, USP (Normal Saline I.V. Flush Syringe) is available in different volumes. It is sterile; contain no preservative, non-pyrogenic aqueous isotonic solution of Sodium Chloride. This DEVICE is not made with natural rubber latex. This DEVICE is not made with DEHP. Each mL contains 9mg/mL of sodium chloride. The total osmolar concentration of solution is 0.308 mOsm/mL (calc.). The pH of the solution is 4.5 – 7.0. The luer lock syringes are made using polypropylene.

Flush Syringe: The flush syringes are steam sterilized prior to final packaging. **DO NOT PLACE THE FLUSH SYRINGE ON STERILE FIELD.**

Sterile Field Ready: The sterile field ready syringes are labeled and packaged in pouch and then steam sterilized. **THESE SYRINGES CAN BE USED IN STERILE FIELD.**

CLINICAL PHARMACOLOGY: 0.9% Sodium Chloride Injection, USP (Normal Saline I.V. Flush Syringe) is a sterile aqueous solution having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to the tissues.

INDICATIONS AND USAGE: For use as a sterile isotonic injection for flushing compatible intravenous tubing and / or indwelling access devices only. The flush syringes are not to be used for reconstitution or dilution of medication.

PRECAUTIONS: Caution must be exercised to avoid any contact between Normal Saline Flush Solution and incompatible drugs. Consult appropriate literature for compatibility data before flushing the catheter containing drug product, with 0.9% Sodium Chloride Injection, USP. Do not administer if a precipitate appears, or if solution is discolored, cloudy, or hazy, or if the syringe is damaged. For single use only. Discard unused portion.

ADVERSE REACTIONS: None known from use in these quantities.


DRUG ABUSE AND DEPENDENCE: Not Applicable.

OVERDOSAGE: None known from use in these quantities.

DOSAGE AND ADMINISTRATION: Lock Flush Solution is injected as a single dose via the injection side of the vascular access device. The volume injected is equivalent to that of the vascular access device. After each time the vascular access device has been used for the injection or infusion of medication, the Lock Flush Solution should be injected to restore the patency of the access device. Aspirate before administering any solution via the access device lock in order to confirm patency and location of the catheter tip.

STORAGE: Do not use if solution is discolored, cloudy, hazy, or contains a precipitate, or if the syringe is damaged. Store at 25°C (77°F); excursions permitted to 15° – 30°C (59° – 86°F). Do not freeze.

WARNING: Re-use of single – use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

For detailed handling information see Directions for Use section.

DIRECTIONS FOR USE:

FLUSH SYRINGE

1. Solution and fluid path are sterile and non-pyrogenic if the tip cap is in place, syringe is intact and there is no evidence of leakage. Use proper aseptic technique.
2. Inspect plastic wrapping. Do not use if plastic wrapping is damaged or not intact. (Figure 1)
3. Remove plastic packaging by tearing along perforation. (Figure 2)
4. Do not use if, solution is discolored, cloudy, hazy, or contains a precipitate, or if the syringe is damaged. (Figure 3)
5. With the tip cap of the syringe on, press the syringe forward to properly activate the syringe. Improper activation may difficult syringe use and/or may cause plunger rod separation. **Never draw back rod because the product may become contaminated.** (Figure 4)
6. Remove tip cap. Hold the syringe unit upright and prime to expel any air bubble if present. (Figure 5)
7. Syringe is now ready to use. Per Institution protocol attach flush syringe to access device and flush. (Figure 6)
8. When attempting to aspirate tubing or indwelling device by pulling back on the syringe plunger, while attached to the tubing or indwelling device, use two-handed technique with one hand on the syringe barrel and the other on the plunger. Pull the plunger straight back. Do not pull or bend the plunger sideways. (Figure 7)
9. Use in accordance with intravenous tubing or indwelling device manufacturer's recommendation. (Figure 8)
10. Discard empty unit after use. Discard any unused portion. Do not reuse disposable syringes. (Figure 9)

STERILE FIELD READY SYRINGE

1. Inspect paper pouch. Do not use if the pouch is damaged, the cap is not completely attached to the syringe or not intact.
2. Open paper pouch by opening along the seal and remove the syringe.
3. Follow steps 4 – 10 as mentioned above under FLUSH SYRINGE.

HOW SUPPLIED: Sodium Chloride Injection, USP is available as follows:

Flush Syringe:

Product No.	Item No.	
MIS-1121	64253-111-21	Normal Saline I.V. Flush Syringe; 1 mL fill in 6 mL syringe, LL
MIS-1122	64253-111-22	Normal Saline I.V. Flush Syringe; 2 mL fill in 6 mL syringe, LL
MIS-1152	64253-111-52	Normal Saline I.V. Flush Syringe; 2.5 mL fill in 6 mL syringe, LL
MIS-1123	64253-111-23	Normal Saline I.V. Flush Syringe; 3 mL fill in 6 mL syringe, LL
MIS-1133	64253-111-33	Normal Saline I.V. Flush Syringe; 3 mL fill in 12 mL syringe, LL
MIS-1125	64253-111-25	Normal Saline I.V. Flush Syringe; 5 mL fill in 6 mL syringe, LL
MIS-1135	64253-111-35	Normal Saline I.V. Flush Syringe; 5 mL fill in 12 mL syringe, LL
MIS-1130	64253-111-30	Normal Saline I.V. Flush Syringe; 10 mL fill in 12 mL syringe, LL

The above products are available in boxes of 30, 60 or 120 counts each.

Sterile Field Ready:

Product No.	Item No.	
MIS-1171	64253-111-71	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 1 mL fill in 6 mL syringe, LL
MIS-1172	64253-111-72	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 2 mL fill in 6 mL syringe, LL
MIS-1177	64253-111-77	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 2.5 mL fill in 6 mL syringe, LL
MIS-1173	64253-111-73	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 3 mL fill in 6 mL syringe, LL
MIS-1175	64253-111-75	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 5 mL fill in 6 mL syringe, LL
MIS-1183	64253-111-83	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 3 mL fill in 12 mL syringe, LL
MIS-1185	64253-111-85	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 5 mL fill in 12 mL syringe, LL
MIS-1180	64253-111-80	Normal Saline I.V. Flush Syringe, Sterile Field Ready, 10 mL fill in 12 mL syringe, LL

Sterile Field Ready product is available in boxes of 50 counts each.

CAUTION: Federal (USA) law restricts this device to sale by or on order of a physician.

Medefil, Inc
0,9% Chlorure de Sodium injectable, USP
"Une solution saline normale I.V. RINCER SERINGUE"

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PHARMACODÉPENDANCE: Non applicable.

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