



MEDEFIL, INC.

Description

Normal Saline I.V. Flush Syringe, 5 mL fill in 6 mL syringe

Key Product Features

Brand	Normal Saline I.V. Flush Syringe
Product Description	Normal Saline I.V. Flush Syringe, 5 mL fill in 6 mL syringe
Catalog Number	MIS-1125
NDC / NHRIC	64253011125
DNR (Latex) free	Yes, Not made with Natural Rubber Latex
Sterilization	Terminally Sterilized
Sterile	Solution and Pathway are Sterile
Tip Cap Color	White
Storage	25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze.
Preservatives	No Preservatives
DEHP	DEHP Free
Concentration	0.9% Sodium Chloride (9 mg/mL)
HPIS Classification	580_160_0_0
UNSPSC Category Code	42142619
GTIN #	1 03 6425311125 0
EAN UPC	64253011125



Shipping and Handling info

Packaging Description	60/BX, 960/CASE
Country of Origin	USA
Shelf Life	24 Months

Product Packaging Info.

Packaging Level	DISPENSER BOX (DB)	MASTER CARTON	PALLET
Quantity	1	16	320
Saleable	FALSE	TRUE	TRUE
Length	11 inch.	24 inch.	48 inch.
Width	7 inch.	15 inch.	40 inch.
Height	5 inch.	20 inch.	84 inch.
Weight	2.5 lbs	41 lbs	860 lbs