



Prefilled Heparin Syringes (1 Unit/mL) for Neonatal Intensive Care Units

- Low-dose prefilled heparin syringes designed specifically for use in Neonatal Intensive Care Units (NICU)
- Tip cap and label are color-coded by concentration to decrease medication errors
- Two year shelf life, compared to pharmacy-compounded alternatives (typical shelf life of seven days)
- No Reflux Technology
- Lowest PSI of any Flush Syringe (12mL or less) creates smooth flow of flush solution with less pressure.
- Thicker wall construction eliminates barrel distortion, keeping fluid in the sterile zone of the syringe.
- Contains no preservatives, not made with natural rubber latex & DEHP
- Bar code & NDC number on each syringe




MEDEFIL, INC.

Phone (630) 682-4600 | www.medefilinc.com

Normal Saline I.V. Flush Syringes

Catalog #	Description	NDC Number	Box / Case Quantity
MIS-1123	3mL fill in 6mL Syringe	64253-111-23	60/Box, 960/case
MIS-1125	5mL fill in 6mL Syringe	64253-111-25	60/Box, 960/case
MIS-1133	3mL fill in 12mL Syringe	64253-111-33	60/Box, 960/case
MIS-1135	5mL fill in 12mL Syringe	64253-111-35	60/Box, 960/case
MIS-1130	10mL fill in 12mL Syringe	64253-111-30	60/Box, 960/case



1 Unit/mL Heparin I.V. Flush Syringes

Catalog #	Description	NDC Number	Box / Case Quantity
MIH-4422	2mL fill in 6mL Syringe	64253-444-22	60/Box, 960/case
MIH-4425	5mL fill in 6mL Syringe	64253-444-25	60/Box, 960/case

Heparin I.V. Flush Syringes in 0.9% Sodium Chloride



10 Units/mL Heparin I.V. Flush Syringes

Catalog #	Description	NDC Number	Box / Case Quantity
MIH-2221	1mL fill in 6mL Syringe	64253-222-21	60/Box, 960/case
MIH-2223	3mL fill in 6mL Syringe	64253-222-23	60/Box, 960/case
MIH-2233	3mL fill in 12mL Syringe	64253-222-33	60/Box, 960/case
MIH-2235	5mL fill in 12mL Syringe	64253-222-35	60/Box, 960/case

Heparin I.V. Flush Syringes in 0.9% Sodium Chloride



100 Units/mL Heparin I.V. Flush Syringes

Catalog #	Description	NDC Number	Box / Case Quantity
MIH-3321	1mL fill in 6mL Syringe	64253-333-21	60/Box, 960/case
MIH-3323	3mL fill in 6mL Syringe	64253-333-23	60/Box, 960/case
MIH-3333	3mL fill in 12mL Syringe	64253-333-33	60/Box, 960/case
MIH-3335	5mL fill in 12mL Syringe	64253-333-35	60/Box, 960/case

Heparin I.V. Flush Syringes in 0.9% Sodium Chloride



MEDEFIL, INC.

250 Windy Point Drive

Glendale Heights Illinois 60139

Phone (630) 682-4600 | Fax (630) 681-9100 | www.medefilinc.com