




MEDEFIL, INC.

Description

Heparin IV Flush Syringe; 300 USP Units/3mL (100 USP Units/mL); 3 mL fill in 6 mL Syringe

Key Product Features

Brand	Heparin I.V. Flush Syringe, USP; 100 Units/mL
Product Description	Heparin IV Flush Syringe; 300 USP Units/3mL (100 USP Units/mL); 3 mL fill in 6 mL Syringe
Catalog Number	MIH-3323
NDC / NHRIC	64253033323
DNR (Latex) free	Yes, Not made with Natural Rubber Latex
Sterilization	Terminally Sterilized
Sterile	Solution and Pathway are Sterile
Tip Cap Color	Yellow
Storage	25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze.
Preservatives	No Preservatives
DEHP	DEHP Free
Concentration	100 Units/mL Heparin Sodium in 9 mg/mL (0.9%) Sodium Chloride
HPIS Classification	580_160_0_0
UNSPSC Category Code	42142619
GTIN #	1 03 6425333323 2
EAN UPC	64253033323



Shipping and Handling info

Packaging Description	60/BX, 960/CASE
Country of Origin	USA
Shelf Life	24 Months

Product Packaging Info.

Packaging Level	DISPENSER BOX (DB)	MASTER CARTON	PALLET
Quantity	1	16	320
Saleable	FALSE	TRUE	TRUE
Length	11 inch.	24 inch.	48 inch.
Width	7 inch.	15 inch.	40 inch.
Height	5 inch.	20 inch.	84 inch.
Weight	2.3 lbs	37 lbs	780 lbs