



# MEDEFIL, INC.

## Description

Heparin I.V. Flush Syringe, USP; 300 Units/3 mL; 3 mL fill in 12 mL Syringe

## Key Product Features

Brand	Heparin I.V. Flush Syringe, USP; 100 Units/mL
Product Description	Heparin I.V. Flush Syringe, USP; 300 Units/3 mL; 3 mL fill in 12 mL Syringe
Catalog Number	MIH-3333
NDC / NHRIC	64253033333
DNR (Latex) free	Yes, Not made with Natural Rubber Latex
Sterilization	Terminally Sterilized
Sterile	Solution and Pathway are Sterile
Tip Cap Color	Yellow
Storage	25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze.
Preservatives	No Preservatives
DEHP	DEHP Free
Concentration	100 Units/mL Heparin Sodium in 9 mg/mL (0.9%) Sodium Chloride
HPIS Classification	580_160_0_0
UNSPSC Category Code	42142619
GTIN #	1 03 6425333333 1
EAN UPC	64253033333



## Shipping and Handling info

Packaging Description	60/BX, 960/CASE
Country of Origin	USA
Shelf Life	24 Months

## Product Packaging Info.

Packaging Level	DISPENSER BOX (DB)	MASTER CARTON	PALLET
Quantity	1	16	320
Saleable	FALSE	TRUE	TRUE
Length	11 inch.	24 inch.	48 inch.
Width	7 inch.	15 inch.	40 inch.
Height	5 inch.	20 inch.	84 inch.
Weight	2.5 lbs	40 lbs	840 lbs