



**MEDEFIL, INC.**

### Description

Heparin Lock Flush Syringe; 30 USP Units/3 mL (10 USP Units/mL); 3 mL fill in 12 mL Syringe

### Key Product Features

Brand	Heparin Lock Flush Syringe, USP; 10 Units/mL
Product Description	Heparin Lock Flush Syringe; 30 USP Units/3 mL (10 USP Units/mL); 3 mL fill in 12 mL Syringe
Catalog Number	MIH-2233
NDC / NHRIC	64253-0222-33
DNR (Latex) free	Yes, Not made with Natural Rubber Latex
Sterilization	Terminally Sterilized
Sterile	Solution and Pathway are Sterile
Tip Cap Color	Blue
Storage	25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze.
Preservatives	No Preservatives
DEHP	DEHP Free
Concentration	10 Units/mL Heparin Sodium in 9 mg/mL (0.9%) Sodium Chloride
HPIS Classification	580_160_0_0
UNSPSC Category Code	42142619
GTIN # (Syringe)	10364253222334
EAN UPC	64253-0222-33



### Shipping and Handling info

Packaging Description	60/BX, 960/CASE
Country of Origin	USA
Shelf Life	24 Months

### Product Packaging Info.

Packaging Level	DISPENSER BOX (DB)	MASTER CARTON	PALLET
GTIN #	30364253222338	50364253222332	
Quantity	1	16	320
Saleable	FALSE	TRUE	TRUE
Length	11 inch.	24 inch.	48 inch.
Width	7 inch.	15 inch.	40 inch.
Height	5 inch.	20 inch.	84 inch.
Weight	2.4 lbs	39 lbs	820 lbs