




MEDEFIL, INC.

Description

Normal Saline IV Flush Syringe [0.9% Sodium Chloride Injection]; 5 mL fill in 12 mL Syringe

Key Product Features

| | |
|----------------------|---------------------------------------------------------------------------------------------|
| Brand | Normal Saline I.V. Flush Syringe |
| Product Description | Normal Saline IV Flush Syringe [0.9% Sodium Chloride Injection]; 5 mL fill in 12 mL Syringe |
| Catalog Number | MIS-1135 |
| NDC / NHRIC | 64253-0111-35 |
| DNR (Latex) free | Yes, Not made with Natural Rubber Latex |
| Sterilization | Terminally Sterilized |
| Sterile | Solution and Pathway are Sterile |
| Tip Cap Color | White |
| Storage | 25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze. |
| Preservatives | No Preservatives |
| DEHP | DEHP Free |
| Concentration | 0.9% Sodium Chloride (9 mg/mL) |
| HPIS Classification | 580_160_0_0 |
| UNSPSC Category Code | 42142619 |
| GTIN # (Syringe) | 10364253111355 |
| EAN UPC | 64253-0111-35 |



Shipping and Handling info

| | |
|-----------------------|-----------------|
| Packaging Description | 60/BX, 960/CASE |
| Country of Origin | USA |
| Shelf Life | 24 Months |

Product Packaging Info.

| Packaging Level | DISPENSER BOX (DB) | MASTER CARTON | PALLET |
|-----------------|--------------------|----------------|----------|
| GTIN # | 30364253111359 | 50364253111353 | |
| Quantity | 1 | 16 | 320 |
| Saleable | FALSE | TRUE | TRUE |
| Length | 11 inch. | 24 inch. | 48 inch. |
| Width | 7 inch. | 15 inch. | 40 inch. |
| Height | 5 inch. | 20 inch. | 84 inch. |
| Weight | 2.8 lbs | 44 lbs | 920 lbs |