



MEDEFIL, INC.

Description

Heparin Lock Flush Syringe; 2 USP Units/2 mL (1 USP Unit/mL); 2 mL fill in 6 mL Syringe

Key Product Features

Brand	Heparin Lock Flush Syringe, USP; 1 Unit/mL
Product Description	Heparin Lock Flush Syringe; 2 USP Units/2 mL (1 USP Unit/mL); 2 mL fill in 6 mL Syringe
Catalog Number	MIH-4422
UDI	64253-0444-22
DNR (Latex) free	Not made with Natural Rubber Latex
Sterilization	Terminally Sterilized
Sterile	Solution and Pathway are Sterile
Tip Cap Color	Green
Storage	25°C (77°F), excursions permitted to 15° - 30°C (59° - 86°F). Do not Freeze.
Preservatives	No Preservatives
DEHP	Not made with DEHP
Concentration	1 Unit/mL Heparin Sodium in 9 mg/mL (0.9%) Sodium Chloride
HPIS Classification	580_160_0_0
UNSPSC Category Code	42142619
GTIN # (Syringe)	10364253444224
EAN UPC	64253-0444-22



Shipping and Handling info

Packaging Description	60/BX, 960/CASE
Country of Origin	USA
Shelf Life	24 Months

Product Packaging Info.

Packaging Level	DISPENSER BOX (DB)	MASTER CARTON	PALLET
GTIN #	30364253444228	50364253444222	
Quantity	1	16	320
Saleable	FALSE	TRUE	TRUE
Length	11 inch.	24 inch.	48 inch.
Width	7 inch.	15 inch.	40 inch.
Height	5 inch.	20 inch.	84 inch.
Weight	2.3 lbs	36 lbs	760 lbs